



A cold fever early know A second temperature measurement Care for the baby

## Non-contact infrared forehead thermometer

**Fast Accurate**

*Rapid temperature Care for the baby*



## Can measure the human body, object, indoor temperature

MORE THAN JUST MEASURING TEMPERATURE



**Water temperature**



**Milk temperature**



**Food temperature**



**Room temperature**

/

# Mercury unexpected harm you

It is very important to choose a good thermometer



**The disadvantage of mercury**

# Measurement accuracy

ACCURATE MEASUREMENT IS THE STANDARD

## Baby thermometer features:

**Screen:** The large screen reading

**Fast:** 1 second rapid temperature

**Security:** Avoid cross infection

**Accurate:** Precise measurement

**Backlight:** Three color backlighting at night



# Three color backlighting

Child/adult body temperature gun

Normal body temperature  
32.0°C~37.5°C

Slight fever  
37.6°C~38.5°C

High body temperature  
38.6°C~42.9°C



# Product features

THE PRODUCT FEATURES



01

**Rapid measurement**

**Accurate error-free**

**A second rapid temperature measurement  
Infrared temperature measurement sensors**



**No need contact**

02

# Chip accurate

Automatic measurement. Accurate measurement





03

## Non-contact measurement

Don't need to contact, clean and health, avoid cross infection and fast measurement to solve the problem baby love moving



04

## Human body engineering

Human body engineering design,  
strengthen his hand grip handle experience,  
and be convenient to measure



# Three color optional

THREE COLOR OPTIONAL



# Product parameters

## PRODUCT PARAMETERS



Range	Body	32°C~42.9°C, 85.9°F~109.2°F
	Object	0°C~100°C, 32°F~212°F
Accuracy	Body	±0.2°C/±0.4°F
	Object	±1.0°C/±1.8°F
Resolution	0.1°C/°F	
Working condition	10°C~35°C (50°F~95°F) RH≤85%	
Storage condition	-20°C~55°C (-4°F~131°F) RH≤93%	
Power supply	D.C. 3V 2pcs AAA Batteries	
Power level indicator	Indication for low power level (<2.5V)	
Memory	Automatic recall of 32 memories	
Display	LCD back-lit display	
Reading scale	Celsius or Fahrenheit	
Automatic shut off	In 18 seconds	
Measuring time	0.5 seconds	

# Product analysis

PRODUCT ANALYSIS



## Product analysis



### Accurate probe

Accurate temperature measuring probe. Measurement accuracy. Little error



### Backlit display

Backlit design, temperature can also clearly see in the night for baby, and grandma and grandpa can see without glasses.



### PP material

Environmental protection PP material lets you use the product is safe and assured.



### One-click button

One-click boot measurement button



# Infrared forehead temperature gun

Care for your family





□□□□ □□

**1. □□□□ □□ □□ □□□□□?**

□□, □□□□ 19 □ □□□ □□□□ □□ □□ □□ □□□□□.

**2. □□ □□□□ □□ □ □ □□□□? □□□□ □□**

□□ □□ □ □□ □□□□□ □□□□ □ □□□ □□ □□□□□□□□. □□ □□ □ □□□ □□□□□□. 5 □ □□□ □□□□ □□ □ □□□□.

**□□ □□□□ □□□□□?**

T / T : 30 % □□, □□ □□ □□□□ 70 % □□

**□□□□□□ □□□□ □□□□□?**

ETD : □□ □□ □ 25-35 □ □□□

**5. □□□□ □□□□ □□□□ □□ □ □ □□□□?**

I. □ □□□ □□□□ □□ □□□□ □□ □□□□ □□□□□□□□. □□□□□□ AQL 2.5 // 4.0 □□

B. □□□□ □□ □□ □□ □□□□ □□□□□□□□. □□□□ □□ □□□□ □□□□□□ □□□□ □□ □□□□□□□□.

III. □□□□ □□□□ □□□□ AQL□ □□□□ □□□□ □□□□ □□□□ □□□□□□□□.



**LOTUS UMBRELLA**

☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐

☐☐ / ☐☐ 0086 13335917698

☐☐ ☐☐

[☐☐☐ ☐☐☐☐ ☐☐ ☐☐ ☐☐☐☐☐☐☐☐☐☐](#)