

A cold fever early know A second temperature measurement Care for the baby

Non-contact infrared forehead thermometer

Fast Accurate

Rapid temperature Care for the baby



Can measure the human body, object, indoor temperature

MORE THAN JUST MEASURING TEMPERATURE



Water temperature



Milk temperature



Food temperature



Room temperature

/

Mercury unexpected harm you

It is very important to choose a good thermometer



The disadvantage of mercury

Measurement accuracy

ACCURATE MEASUREMENT IS THE STANDARD

Baby thermometer features:

Screen: The large screen reading

Fast: 1 second rapid temperature

Security: Avoid cross infection

Accurate: Precise measurement

Backlight: Three color backlighting at night



Three color backlighting

Child/adult body temperature gun

Normal body temperature
32.0°C~37.5°C

Slight fever
37.6°C~38.5°C

High body temperature
38.6°C~42.9°C



Product features

THE PRODUCT FEATURES



01

Rapid measurement

Accurate error-free

**A second rapid temperature measurement
Infrared temperature measurement sensors**



No need contact

02

Chip accurate

Automatic measurement. Accurate measurement



03

Non-contact measurement

Don't need to contact, clean and health, avoid cross infection and fast measurement to solve the problem baby love moving



04

Human body engineering

Human body engineering design,
strengthen his hand grip handle experience,
and be convenient to measure



Three color optional

THREE COLOR OPTIONAL



Product parameters

PRODUCT PARAMETERS



Range	Body	32°C~42.9°C, 85.9°F~109.2°F
	Object	0°C~100°C, 32°F~212°F
Accuracy	Body	±0.2°C/±0.4°F
	Object	±1.0°C/±1.8°F
Resolution		0.1°C/°F
Working condition		10°C~35°C (50°F~95°F) RH≤85%
Storage condition		-20°C~55°C (-4°F~131°F) RH≤93%
Power supply		D.C. 3V 2pcs AAA Batteries
Power level indicator		Indication for low power level (<2.5V)
Memory		Automatic recall of 32 memories
Display		LCD back-lit display
Reading scale		Celsius or Fahrenheit
Automatic shut off		In 18 seconds
Measuring time		0.5 seconds

Product analysis

PRODUCT ANALYSIS



Product analysis



Accurate probe

Accurate temperature measuring probe. Measurement accuracy. Little error



Backlit display

Backlit design, temperature can also clearly see in the night for baby, and grandma and grandpa can see without glasses.



PP material

Environmental protection PP material lets you use the product is safe and assured.



One-click button

One-click boot measurement button



Infrared forehead temperature gun

Care for your family





□□□□ □□

1. □□□□ □□ □□ □□□□□?

□□, □□□□ 19 □ □□□ □□□□ □□ □□ □□ □□□□□.

2. □□ □□□□ □□ □ □ □□□□? □□□□ □□

□□ □□ □ □□ □□□□□ □□□□ □ □□□ □□ □□□□ □□□□. □□ □□ □ □□□ □□□□□. 5 □ □□□ □□□□ □□ □ □□□□.

□□ □□□□ □□□□□?

T / T : 30 % □□, □□ □□ □□□□ 70 % □□

□□□□□□ □□□□ □□□□□?

ETD : □□ □□ □ 25-35 □ □□□□

5. □□□□ □□□□ □□□□ □□ □ □ □□□□?

I. □ □□□ □□□□ □□ □□□□ □□ □□□□ □□□□□. □□□□□□ AQL 2.5 // 4.0 □□

B. □□□□ □□ □□ □□ □□□□ □□□□□. □□□□ □□ □□□□ □□□□□ □□□□ □□ □□□□□.

III. □□□□ □□□□ □□□□ AQL□ □□□□ □□□□ □□□□ □□□□ □□□□□.



LOTUS UMBRELLA

☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐

☐☐ / ☐☐ 0086 13335917698

☐☐ ☐☐

[☐☐☐ ☐☐☐☐ ☐☐ ☐☐ ☐☐☐☐☐☐☐☐☐☐](#)